

XPERT HAND.



SET DESCRIPTION



XPert HAND.

Intended purpose:

The implants of the Xpert Hand range are intended for fixation of hand fractures, osteotomies and arthrodeses in adults.

Contraindications:

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.

3 CONFIGURATIONS FOR THE XPert HAND SET


1 • COMPLETE SET



- General instrumentation (including optional pliers)
- Phalangeal module (instrumentation, plates & screws)
- Metacarpal module (instrumentation, plates & screws)

Pages 3 & 4

2 • PHALANGEAL SET (2 LEVELS)



- General instrumentation (including optional pliers)
- Phalangeal module (instrumentation, plates & screws)

Page 5

3 • METACARPAL SET (2 LEVELS)




- General instrumentation (including optional pliers)
- Metacarpal module (instrumentation, plates & screws)

Page 6

1. Complete set (1ST LEVEL)



B3: General instrumentation
(ANC1868/B3)




GENERAL INSTRUMENTATION FOR THE FULL RANGE

ANC1613
Bending pliers (x2) -
For 0.8 & 1.3 mm plates



ANC1612
Cutting pliers -
For 0.8 & 1.3 mm plates



OPTIONAL RANGE OF FORCEPS AND PLIERS

ANC1609
Retractor forceps -
For 0.8 & 1.3 mm plates



ANC1611
Plate holding forceps -
For 0.8 & 1.3 mm plates




ANC503
150 mm reduction
forceps



ANC504
150 mm pointed
reduction forceps





1. Complete set (2ND LEVEL)



M2: Metacarpal module

(ANC1868/M2)

Instrumentation, 1.3 mm plates and screws for metacarpals
(see page 10 & 11 for the references)





1. Complete set (3RD LEVEL)



M1: Phalangeal module


(ANC1868/M1)

Instrumentation, 0.8 mm plates and screws for phalanges
(see page 10 & 11 for the references)




2. Phalangeal set (1ST LEVEL)

B3: General instrumentation (ANC1868/B3)




GENERAL INSTRUMENTATION FOR THE FULL RANGE

ANC1613
Bending pliers (x2) -
For 0.8 & 1.3 mm plates



ANC1612
Cutting pliers -
For 0.8 & 1.3 mm plates



OPTIONAL RANGE OF FORCEPS AND PLIERS

ANC1609
Retractor forceps -
For 0.8 & 1.3 mm plates



ANC1611
Plate holding forceps -
For 0.8 & 1.3 mm plates



ANC503
150 mm reduction
forceps



ANC504
150 mm pointed
reduction forceps





2. Phalangeal set (2ND LEVEL)

M1: Phalangeal module


(ANC1868/M1)

Instrumentation, 0.8 mm plates and screws for phalanges
(see page 10 & 11 for the references)




3. Metacarpal set (1ST LEVEL)

B3: General instrumentation (ANC1868/B3)




GENERAL INSTRUMENTATION FOR THE FULL RANGE

ANC1613
Bending pliers (x2) -
For 0.8 & 1.3 mm plates



ANC1612
Cutting pliers -
For 0.8 & 1.3 mm plates



OPTIONAL RANGE OF FORCEPS AND PLIERS

ANC1609
Retractor forceps -
For 0.8 & 1.3 mm plates



ANC1611
Plate holding forceps -
For 0.8 & 1.3 mm plates



ANC503
150 mm reduction
forceps





ANC504
150 mm pointed
reduction forceps



3. Metacarpal set (2ND LEVEL)

M2: Metacarpal module (ANC1868/M2)

Instrumentation, 1.3 mm plates and screws for metacarpals
(see page 10 & 11 for the references)



Implant references not presented above.

Ø1.2 mm locking screws*

Ref.	Description
SAT1.2L06D to SAT1.2L16D	Ø1.2 mm locking screw - L06 to 16 mm (1mm increments)
SAT1.2L18D to SAT1.2L20D	Ø1.2 mm locking screw - L18 to 20 mm (2mm increments)



*Non anodized

Ø1.2 mm non-locking screws*

Ref.	Description
CAT1.2L06D to CAT1.2L16D	Ø1.2 mm non-locking screw - L06 to 16 mm (1mm increments)
CAT1.2L18D to CAT1.2L20D	Ø1.2 mm non-locking screw - L18 to 20 mm (2mm increments)



*Light pink anodized

Ø1.5 mm locking screws*

Ref.	Description
SAT1.5L06D to SAT1.5L16D	Ø1.5 mm locking screw - L06 to 16 mm (1mm increments)
SAT1.5L18D to SAT1.5L24D	Ø1.5 mm locking screw - L18 to 24 mm (2mm increments)



*Non anodized

Ø1.5 mm non-locking screws*


Ref.	Description
CAT1.5L06D to CAT1.5L16D	Ø1.5 mm non-locking screw - L06 to 16 mm (1mm increments)
CAT1.5L18D to CAT1.5L24D	Ø1.5 mm non-locking screw - L18 to 24 mm (2mm increments)



*Light pink anodized

Ø2.0 mm locking screws*


Ref.	Description
SAT2.0L06D to SAT2.0L16D	Ø2.0 mm locking screw - L06 to 16 mm (1mm increments)
SAT2.0L18D to SAT2.0L24D	Ø2.0 mm locking screw - L18 to 24 mm (2mm increments)



*Non anodized

Ø2.0 mm non-locking screws*

Ref.	Description
CAT2.0L06D to CAT2.0L16D	Ø2.0 mm non-locking screw - L06 to 16 mm (1mm increments)
CAT2.0L18D to CAT2.0L24D	Ø2.0 mm non-locking screw - L18 to 24 mm (2mm increments)



*Light pink anodized

Ø2.3 mm locking screws*


Ref.	Description
SAT2.3L06D to SAT2.3L16D	Ø2.3 mm locking screw - L06 to 16 mm (1mm increments)
SAT2.3L18D to SAT2.3L24D	Ø2.3 mm locking screw - L18 to 24 mm (2mm increments)



*Non anodized

Ø2.3 mm non-locking screws*

Ref.	Description
CAT2.3L06D to CAT2.3L16D	Ø2.3 mm non-locking screw - L06 to 16 mm (1mm increments)
CAT2.3L18D to CAT2.3L24D	Ø2.3 mm non-locking screw - L18 to 24 mm (2mm increments)




*Light pink anodized

Implant references not presented above.

0.8 mm plates

Ref.	Description	Bendable	Cuttable
PTSM12D	0.8 mm hybrid plate - 12 holes	Yes	Yes
PTSS6D	0.8 mm straight plate - 6 holes	Yes	Yes
PTSS10D	0.8 mm straight plate - 10 holes	Yes	Yes
PTSD8D	0.8 mm glide plate - 8 holes	Yes	No
PTSD10D	0.8 mm glide plate - 10 holes	Yes	No
PTDP6D	0.8 mm phalangeal head plate - Medial for right hand / Lateral for left hand - 6 holes	Yes	No
PTGP6D	0.8 mm phalangeal head plate - Medial for left hand / Lateral for right hand - 6 holes	Yes	No
PTSH1D	0.8 mm hook plate - 1 hole	No	No



1.3 mm plates

Ref.	Description	Bendable	Cuttable
MTSM12D	1.3 mm hybrid plate - 12 holes	Yes	Yes
MTSS6D	1.3 mm straight plate - 6 holes	Yes	Yes
MTSS10D	1.3 mm straight plate - 10 holes	Yes	Yes
MTSD8D	1.3 mm glide plate - 8 holes	Yes	No
MTSD10D	1.3 mm glide plate - 10 holes	Yes	No
MTSH6D	1.3 mm Rolando hook plate - 6 holes	Yes	No
MTSR7D	1.3 mm rotation plate - 7 holes	Yes	No
MTSO6D	1.3 mm double oblong straight plate - 6 holes	Yes	No



Trays references.


Complete set

Ref.	Description
ANC1868/C	Xpert Hand set - Lid
ANC1868/M2	Xpert Hand set - Module 2 - 1.3 mm plates
ANC1868/M1	Xpert Hand set - Module 1 - 0.8 mm plates
ANC1868/B3	Xpert Hand set - Base - 3 levels



Phalangeal set

Ref.	Description
ANC1868/C	Xpert Hand set - Lid
ANC1868/M1	Xpert Hand set - Module 1 - 0.8 mm plates
ANC1868/B2	Xpert Hand set - Base - 2 levels



Metacarpal set

Ref.	Description
ANC1868/C	Xpert Hand set - Lid
ANC1868/M2	Xpert Hand set - Module 2 - 1.3 mm plates
ANC1868/B2	Xpert Hand set - Base - 2 levels



This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions including cleaning and sterilization before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

Manufacturer: Newclip Technics - Compo Kit EN - Xpert Hand - Ed.1 - 07/2026 - Medical devices: class IIb - CE1639 SGS BE - US Class: II - Read labelling and instructions before the use of Newclip Technics medical devices. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. Non-contractual pictures.
Newclip Technics - 45 rue des Garottières - 44115 Haute Goulaine, France. Our subsidiaries: Newclip USA - Newclip Australia - Newclip GMBH - Newclip Japan - Newclip Iberia - Newclip Belgium - Newclip Italia - Newclip UK.

newcliptechnics.com

